

### Boost Productivity and Reliability

- Measures and compresses pulses in seconds
- Uses the MIIPS® auto-calibration and auto-compression technology
- Complete pulse shaping solution (includes computer, spectrometer and nonlinear optical detection)
- Most affordable pulse shaping solution
- Small footprint (6-inch cube)
- Winner of the 2009 PhAST / Laser Focus World Innovation Award



### Eliminate Manual Tweaking

*With Push-Button Pulse Characterization*

**Includes the MIIPS® technology:**

MIIPS® is an automated procedure for measurement and compression of optical pulses. It uses a calibrated pulse shaper to introduce a set of reference phase functions and monitors their effect on spectrally resolved nonlinear response such as second harmonic generation.

Mathematical analysis of the recorded spectra provides a direct measurement of high-order pulse dispersion. The measured spectral phase can be compensated by the pulse shaper to compress the laser pulses to their transform limit at the target, without manual tweaking.



## System Specifications


Number of pixels	128
Operating wavelength range*	700 - 900 nm
Maximum spectral window	200 nm
Shaping of spectral phase independent of repetition rate	

\*Contact us if a different wavelength range or spectral window is desired.  
We have delivered systems working in the 500-1800nm range.

## Schematic

Recommended beam diameter ~2 mm


Dimensions L x W x H 150 x 150 x 150 mm (6.0 x 6.0 x 6.0 in.)



Biophotonic Solutions continuously follows a strict product improvement and evaluation program. Specifications are subject to change without notice.

### Example showing accurate shaping before regenerative amplification

Reduces pulse duration from ~50s to ~35fs in most commercial regenerative amplifiers by removing third and fourth order dispersion. Call us to request system specific information.



Copyright (c) 2010 Biophotonic Solutions Inc.